

HPLC Troubleshooting Mini Guide

Peak Issues



Contents

1. Overview
2. Peak Issues
3. Key Problem Areas and Preventive Maintenance
4. Protect Your HPLC Column



Overview

1

Locating and Correcting the Problem

A systematic approach to identifying the problem is the best path to troubleshooting your HPLC system. This guide is organized by four major categories of symptoms to help you quickly identify the source of the problem(s) you are encountering:

- pressure abnormalities
- leaks
- peak problems
- baseline issues



When you have corrected the problem, record the incident in the system recordbook to help with future problems.

Prevention

Many LC problems can be prevented with routine preventive maintenance such as replacing pump seals regularly. Consistent preventive maintenance practices will enhance lab productivity, avoid system critical damage, equipment downtime and costly repairs.

Where to Get Additional Help


1. Chat with Phenomenex technical experts. Phenomenex has experienced technical consultants who can assist you with any chromatography issue in real time. To chat now go to www.phenomenex.com/chat.
2. The operator's and service manuals for the instrument should be consulted. These contain exploded diagrams, troubleshooting procedures for specific models, and part numbers to help you order replacement parts.
3. Other people in the lab may have had experience solving a problem which is giving you trouble; they can be a helpful resource.
4. The manufacturer of your instrument can help you. Most LC manufacturers offer free technical support to their customers.
5. Phenomenex offers seminars on HPLC/UHPLC. Join PhenoAcademy for specialty troubleshooting webinars, www.phenomenex.com/phenoacademy
6. Other resources:
 - J.W. Dolan and L.R. Snyder, **Troubleshooting LC Systems**, Humana Press, NJ (1989).
 - L.R. Snyder and J.J. Kirkland, **Introduction to Modern Liquid Chromatography**, 2nd ed., Wiley, NY (1979).
 - D.J. Runser, **Maintaining and Troubleshooting HPLC Systems - A User's Guide**, Wiley, NY (1981).
 - J.W. Dolan, "LC Troubleshooting", LC/GC Magazine. This is a monthly column.

Peak Issues

2

Many problems in the LC system show up as changes in the chromatogram. Some of these can be solved by changes in the equipment; however, others require modification of the assay procedure. Selecting the proper column type and mobile phase are keys to “good chromatography.”

Peak tailing

POSSIBLE CAUSE	SOLUTION
1. Blocked guard/frit or column contamination	1. a. Replace/remove guard b. Reverse flush column (if allowed) c. Replace column
2. Column void	2. Replace column; avoid sudden pressure shocks. Ensure operating pressure is within column advised operating conditions and that mobile phase pH is within operating guidelines for the column utilized.
3. Interfering peak	3. a. Use longer column b. Change mobile phase and/or column/selectivity
4. Wrong mobile phase pH	4. Adjust pH. For basic compounds, lower pH usually provides more symmetrical peaks
5. Inadequate buffer concentration	5. Increase buffer concentration
6. Sample reacting with active sites	6. a. Use lower pH to deactivate residual silanols b. Use columns specialized in reducing secondary interactions (consult Phenomenex) c. Utilize ion-pair reagent or volatile basic modifier
	
7. Extra-column effect	7. a. Replumb system (shorter, narrower tubing) b. Use smaller volume detector cell c. Check for correct column connection

Peak fronting


POSSIBLE CAUSE	SOLUTION
1. Sample solvent stronger than mobile phase	1. Use mobile phase for injection solvent
2. Sample overload	2. a. Decrease sample concentration or volume b. Increase column ID

Distortion of larger peaks

POSSIBLE CAUSE	SOLUTION
1. Sample overload	1. a. Decrease sample concentration or volume b. Increase column ID

Peak Issues (continued)

All peaks splitting

POSSIBLE CAUSE	SOLUTION
1. Contamination on guard or analytical column inlet 	1. Remove guard column or cartridge and attempt analysis. Replace guard if necessary. If the analytical column is obstructed, reverse and flush. If problem persists, column may be fouled with strongly retained contaminants. Use appropriate restoration procedure. If problem persists, inlet is probably plugged. Replace column.
2. Sample solvent incompatible with mobile phase	2. Change solvent. Whenever possible, inject samples in starting mobile phase.

Some peaks splitting

POSSIBLE CAUSE	SOLUTION
1. Sample solvent incompatible with mobile phase	1. Early eluting peaks may be more affected by sample solvent mismatch. Inject sample in starting mobile phase.
2. Split peak may be two coeluting peaks	2. Use a smaller injection volume to see if peaks become more distinct. Optimize method to increase resolution.
3. pH of the mobile phase is too close to the pKa of the functional group on the compound, or mobile phase is inadequately buffered	3. Adjust and buffer the mobile phase pH +/- 2 units above or below the pKa of the compound's ionizable functional group, and buffer at that pH with an appropriate.

Distortion of early peaks

POSSIBLE CAUSE	SOLUTION
1. Wrong injection solvent	1. a. Reduce injection volume b. Use weaker injection solvent

Extra peaks

POSSIBLE CAUSE	SOLUTION
1. Late-eluting peak from previous injection	1. a. Increase run time or gradient slope b. Increase flow rate
2. Negative or ghost peaks	2. a. Check purity of mobile phase. Mobile phase contaminants commonly accumulate during weak portion of a gradient and elute as strong solvent increases. b. Use mobile phase as injection solvent c. Reduce injection volume
3. Sample Contamination	3. a. Perform additional sample preparation b. Use and/or replace guard column/cartridge c. Clean column

Peak Issues (continued)

Retention time drifts

POSSIBLE CAUSE	SOLUTION
1. Poor temperature control	1. Thermostat column
2. Mobile phase changing	2. a. Prepare fresh mobile phase. Preventative measures include using SecurityCap to prevent solvent evaporation and using an effective buffer concentration. Buffer range capacity is +/- 1 pH from buffer pKa b. Ensure that mobile phase is homogenous
3. Poor column equilibration	3. Allow more time for column equilibration between runs
4. Injection Port Contamination	4. a. Clean Injection Port & flush sample loop b. Replace Needle Wash; increase needle wash strength

Abrupt retention time changes

POSSIBLE CAUSE	SOLUTION
1. Flow rate change	1. Reset flow rate, confirm pump is delivering correct flow
2. Air bubble in pump	2. Bleed air from pump
3. Improper mobile phase	3. a. Replace with proper mobile phase b. Set proper mobile phase mixture on controller
3. Degradation of stationary phase	3. Replace column; confirm pH limits of column

Broad Peaks

POSSIBLE CAUSE	SOLUTION
1. Excessive extra column volume	1. Reduce extra column volume (tubing length/ID, sample loop size, flow cell size, etc.)
2. Sample solvent incompatible with mobile phase	2. Dilute in mobile phase or weak solvent
3. Mobile phase pH incorrect	3. Use buffered mobile phase, commonly +/- 2 pH units from analyte pKa is ideal
4. Incorrect or insufficient buffer	4. Check that the correct buffer is being used. Buffer capacity is +/- 1 pH unit from buffer pKa. Increase buffer concentration.
5. High longitudinal diffusion	5. Retention time too long. Use gradient elution, faster flow rates, stronger eluent, and/or less retentive stationary phase (e.g. core-shell)
6. Column overload	6. a. Reduce sample concentration or volume. b. Use larger column ID
7. Detector response time too long	7. Increase detector scan rate

Key Problem Areas and Preventive Maintenance



The chart below lists the most common problems that occur with each LC module. In the right-hand column are listed preventive maintenance practices that can reduce the failure rate. The numbers in parentheses are suggested intervals between maintenance. The operator's and service manuals for your LC may have additional suggestions for preventive maintenance of your model of LC.

Reservoir

POSSIBLE CAUSE	SOLUTION
1. Blocked inlet frit	1. a. Replace (3-6 mo.) b. Filter mobile phase, 0.5 μ m filter
2. Gas bubbles	2. Degas mobile phase

Pump

POSSIBLE CAUSE	SOLUTION
1. Air bubbles	1. Degas mobile phase
2. Pump seal failure	2. Replace (3 mo.)
3. Check valve failure	3. Filter mobile phase, use inlet-line frit. Keep spare.

Injector

POSSIBLE CAUSE	SOLUTION
1. Rotor seal wear	1. a. Don't overtighten b. Filter samples

Column

POSSIBLE CAUSE	SOLUTION
1. Blocked frit	1. a. Filter mobile phase b. Filter samples c. Use in-line filter and/or guard column
2. Void at head of column	2. a. Avoid mobile phase pH at or near column max pH limit b. Install column at low flow rates to avoid pressure shock c. Use guard column

Protect Your HPLC Column

4


Contaminants Can Cause the Following:

- High Backpressure
- Split Peaks
- Broad Peaks
- Baseline Noise
- Baseline Drift
- Loss of Resolution
- Irreversible Column Damage
- System Damage

Protect Your HPLC Column. Protect Your Results.

The SecurityGuard™ and SecurityGuard ULTRA cartridge systems effectively protect analytical columns from the damaging effect of contaminants that could impact results and data quality. Either cartridge system is designed to trap contaminants without altering your chromatography.

How It Works:



SecurityGuard and SecurityGuard ULTRA standard can adjust to fit any manufacturer's female/inverted endfitting.



Additional information can be found at www.phenomenex.com/securityguard

Phenex™ Syringe Filters

For Sample and Solvent Filtration Prior to Chromatography

- Less system downtime
- More consistent, reproducible results
- Increased column lifetime

Features and Benefits

- Low protein adsorption
- Broad chemical compatibility
- Minimized extractables
- Excellent flow rate
- High total throughput
- Low hold-up volume
- Certified quality
- 100 % integrity tested



MEMBRANE TYPES

RC (Regenerated Cellulose)

PTFE, Teflon® (Polytetrafluoroethylene)

PES (Polyethersulfone)

PVDF (Polyvinylidene Fluoride)

NY (Nylon)

CA (Cellulose Acetate)

GF (Glass Fiber)



To order and find out more about this product visit www.phenomenex.com/phenex

HPLC Troubleshooting Mini Guide

Peak Issues

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.
www.phenomenex.com/behappy

Australia
t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria
t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada
t: +1 (800) 543-3681
info@phenomenex.com

China
t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic
t: +420 272 017 077
cz-info@phenomenex.com

Denmark
t: +45 4824 8048
nordicinfo@phenomenex.com

Finland
t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France
t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany
t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong
t: +852 6012 8162
hkinfo@phenomenex.com

India
t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia
t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland
t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy
t: +39 051 6327511
italiainfo@phenomenex.com

Japan
t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg
t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico
t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands
t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand
t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway
t: +47 810 02 005
nordicinfo@phenomenex.com

Poland
t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal
t: +351 221 450 488
ptinfo@phenomenex.com

Singapore
t: 800-852-3944
sginfo@phenomenex.com

Slovakia
t: +420 272 017 077
sk-info@phenomenex.com

Spain
t: +34 91-413-8613
espinfo@phenomenex.com

Sweden
t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland
t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan
t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand
t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom
t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA
t: +1 (310) 212-0555
info@phenomenex.com

All other countries/regions
Corporate Office USA
t: +1 (310) 212-0555
www.phenomenex.com/chat



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

No part of this booklet may be copied without prior written permission from Phenomenex, Inc. USA. While every attempt has been made to ensure the accuracy of the information contained in this guide, Phenomenex assumes no responsibility for its use. We welcome any additions or corrections for incorporation into future editions. Phenex, SecurityGuard and BE-HAPPY are trademarks of Phenomenex. Teflon is a registered trademark of E.I. du Pont de Nemours and Co. SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362. CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

© 2024 Phenomenex, Inc. All rights reserved.