

Ibuprofen and Paracetamol Tablets – Assay and Related Substances Method as per IP Monograph

Dr. Trivikram Reddy Gundala¹, Swetha Kotikalaupdi¹, and Dr. Bryan Tackett²

¹India PhenoLogix Lab, Mahadeva Kodigehalli, Hobli, Jala Taluka, Bengaluru 562149, India

²Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501-1430, USA

Overview

The combination of Paracetamol and Ibuprofen has been found to be efficacious in a variety of acute pain states, including postoperative pain, dysmenorrhea and musculoskeletal pain. The Paracetamol / Ibuprofen combination is known to provide improved pain relief compared to either drug alone and is usually offered as a substitute to codeine-based analgesics for short-term management of pain.

This application note highlights the use of the Phenomenex Luna® C8(2) column for the identification of Ibuprofen and Paracetamol tablets as per Indian Pharmacopoeia (IP). The requirements to meet system suitability for the assay is a percent relative standard deviation (%RSD) no more than (NMT) 2 % for each peak, and for related substances a resolution of no less than (NLT) 10 and a %RSD NMT 2 % for each of the principal peaks. The method meets the system suitability requirements as per the monograph when analyzing this drug product sample and demonstrates excellent precision over multiple injections.



LC-UV Conditions - Assay Method

Column: Luna® 5 μm C8(2)
Dimension: 250 x 4.6 mm
Part No.: 00G-4249-E0

Mobile Phase: Buffer/Acetonitrile/Methanol (75:24.5:0.5, v/v/v)

Buffer: 3.9 g of Sodium dihydrogen phosphate dihydrate (25 mM) and 8.9 g of Disodium hydrogen phosphate

dihydrate (50 mM) are dissolved in 1000 mL water and adjust pH to 7.0 with Orthophosphoric acid

Flow Rate: 1.5 mL/min (Isocratic)

Injection Volume: 5 μL Temperature: 35 °C

Detection: UV @ 220 nm


LC System: Waters® ACQUITY® H-Class HPLC

Sample Preparation

Transfer powdered sample (approximate weight 0.384 g) of Paracetamol and Ibuprofen tablet which contains 0.16 g Ibuprofen in 100 mL volumetric flask, dissolve and dilute with the diluent (mobile phase) and filter the solution using CLARIFY-NY 13 mm Syringe Filters (AF8-7707-12). Take 10 mL of the above solution in 100 mL volumetric flask and dilute with the diluent.

NOTE: Average tablet weight: 0.963 g / tab Strength of Ibuprofen: 400 mg Strength of Paracetamol: 325 mg

Figure 1. Overlay Chromatogram of Six Injections of Paracetamol and Ibuprofen Test Solution for Assay on Luna 5 µm C8(2)



Compound (Precision, N=6)	Retention Time (min)	Retention Time %RSD	Area	Area %RSD	Symmetry Factor	Resolution
Paracetamol	2.53	0.08	767341	0.18	1.02	-
Ibuprofen	16.91	0.09	1381303	0.11	1.13	48.4

Sample Preparation

Transfer powdered sample (approximate weight 0.241 g) of Paracetamol and Ibuprofen tablet which contains 0.1 g Ibuprofen in 50 mL volumetric flask, add 30 mL of diluent (see below) sonicate for 10 minutes. Cool the solution, add diluent up to volume with diluent, and filter the solution CLARIFY-NY 13 mm Syringe Filters (AF8-7707-12).

Diluent: Solution A/Methanol (40:60, v/v)

Solution A: $6.8\,\mathrm{g}$ of Monobasic potassium phosphate in $1000\,\mathrm{mL}$ water and pH adjusted to $3.0\,\mathrm{with}$ Orthophosphoric acid.

LC-UV Conditions – Related Substances Method

Column: Luna® 5 µm C8(2)

Dimension: 250 x 4.6 mm

Part No.: 00G-4249-E0

Mobile Phase: A: 6.8 g of Monobasic potassium phosphate in 1000

mL water (pH adjusted to 3.0 with Orthophosphoric acid) with 0.25 g Hexane 1-sulphonic acid sodium salt

B: Acetonitrile

	2171001011111	
Gradient:	Time (min)	%В
	0	2
	6	2
	16	20
	22	35
	55	40
	63	55
	63.1	2
	70	2


Flow Rate: 1.5 mL/min Injection Volume: 10 µL

Temperature: 25 °C

Detection: UV @ 214 nm

LC System: Waters® ACQUITY® H-Class HPLC

Figure 2. Paracetamol and Ibuprofen Test Solution for Related Substances on Luna 5 μm C8(2)



Compound	Retention Time (min)	Area	Theoretical Plates	Symmetry Factor	Resolution
Paracetamol	12.60	17747139	106396	1.12	-
Ibuprofen	60.65	29994847	267055	1.09	154.02

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

t: +61 (0)2-9428-6444 auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301 anfrage@phenomenex.com

Belaium

t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681 info@phenomenex.com

t: +86 400-606-8099 cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077 cz-info@phenomenex.com

Denmark

t: +45 4824 8048 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

t: +49 (0)6021-58830-0 anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162 hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400 indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707 indoinfo@phenomenex.com

t: +353 (0)1 247 5405 eireinfo@phenomenex.com

Italy t: +39 051 6327511 italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262 jpinfo@phenomenex.com

Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226 tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700 nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951 nzinfo@phenomenex.com

Norway t: +47 810 02 005 nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72 pl-info@phenomenex.com

Portugal t: +351 221 450 488 ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944 sginfo@phenomenex.com

Slovakia t: +420 272 017 077 sk-info@phenomenex.com

Spain

t: +34 91-413-8613 espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287 thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367 ukinfo@phenomenex.com

t: +1 (310) 212-0555 info@phenomenex.com

All other countries/regions

Corporate Office USA t: +1 (310) 212-0555 info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna is a registered trademark and BE-HAPPY is a trademark of Phenomenex.

Waters and ACQUITY are registered trademarks of Waters Technologies Corporation. Disclaimer

Phenomenex is in no way affiliated with Waters Technologies Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures. © 2021 Phenomenex, Inc. All rights reserved.

