

Performance of the Agela Venusil® 5 μm XBP C18(L) Column

Dawn Chen, PhD and Bryan Tackett, PhD Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

The analysis of an ABN3 sample was used to assess the performance of the Agela Venusil 5 µm XBP C18(L) column compared to the Avantor® ACE® 5 μm C18 column and Avantor ACE Excel® 5 μm C18 column. ABN3 was analyzed because it allows for probing of the surface of the matrix for overall retention and peak shape. This is because the basic analyte is sensitive to the active silanol groups on the matrix surface. Due to the higher column efficiency of the Venusil 5 μm XBP C18(L) column, there was greater resolution and separation of the analytes. There was also greater intensity compared to the ACE and ACE Excel 5 µm C18 columns. In this application note, we show that the Venusil 5 µm XBP C18(L) column outperforms the ACE 5 µm C18 column as well as the ACE Excel 5 µm C18 column under these conditions.

LC Conditions

Column: Venusil 5 μm XBP C18(L) (VX950505-L)

ACE 5 µm C18 ACE Excel 5 µm C18

Dimensions: 50 x 4.6 mm

Mobile Phase: A: 0.5 % Formic acid in water

B: Acetonitrile

Gradient: Time (min) %B 0.25 10 0.5 20 3.5 40 90 3.5 90 10 10

Flow Rate: 1 mL/min Injection Volume: 10 µL Temperature: 20 °C

LC System: Agilent® 1260 Infinity

Detection: MS/MS

Detector: SCIEX® API 4000 QTRAP®

MS Conditions

Ion Source: Positive or Negative

Source Temperature: 600 °C

GS1: 30 **GS2**: 30 **CUR:** 10

IS: +4500, or -4500

Table 1. MRM Transitions

Peak No.	Analyte	Q1 (m/z)	Q3 (m/z)
1	3-Methyl, 4-Nitrobenzoic acid	181.15	137.5
2	Nortriptyline	266	92.5
3	2-Hydroxy, 5-Methylbenzaldehyde	135.7	108.3

Table 2. Column Specifications

Specifications	Venusil 5 μm XBP C18(L)	ACE 5 μm C18	ACE Excel 5 μm C18
Pore Size (Å)	150	100	100
Carbon Load (%)	15	15.5	15.5
Surface Area (m²/g)	200	300	300
Endcapping	Yes	Yes	Yes
pH Range	1.5 – 9.0	2.0 - 8.0	2.0 - 8.0

Figure 1. Separation of ABN3 on the Venusil® 5 μm C18(L) Column




Figure 2. Separation of ABN3 on the ACE® 5 μm C18 Column





Figure 3. Separation of ABN3 on the ACE® Excel® 5 μm C18 Column



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

t: +61 (0)2-9428-6444 auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301 anfrage@phenomenex.com

Belaium

t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681 info@phenomenex.com

t: +86 400-606-8099 cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077 cz-info@phenomenex.com

Denmark

t: +45 4824 8048 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

t: +49 (0)6021-58830-0 anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162 hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400 indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707 indoinfo@phenomenex.com

t: +353 (0)1 247 5405 eireinfo@phenomenex.com

Italy t: +39 051 6327511 italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262 jpinfo@phenomenex.com

Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226 tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700 nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951 nzinfo@phenomenex.com

Norway t: +47 810 02 005 nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72 pl-info@phenomenex.com

Portugal t: +351 221 450 488 ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944 sginfo@phenomenex.com

Slovakia t: +420 272 017 077 sk-info@phenomenex.com

Spain

t: +34 91-413-8613 espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287 thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367 ukinfo@phenomenex.com

t: +1 (310) 212-0555 info@phenomenex.com

All other countries/regions

Corporate Office USA t: +1 (310) 212-0555 info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Avantor, ACE, and Excel are registered trademarks of VWR International, LLC. SCIEX and QTRAP are registered trademarks of AB SCIEX Pte. Ltd. Agilent is a registered trademark of Agilent Technologies, Inc. BE-HAPPY is a trademark of Phenomenex. Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of Phenomenex, Inc. or their respective owners in the United States and/or certain other countries.

Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with VWR International, LLC or Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures. © 2022 Phenomenex, Inc. All rights reserved.





