


Robust Analysis of Nitrosamine Using Strata™ Activated Carbon Extraction

Eduardo Morales¹, Agustin Pierri, PhD¹, Ramkumar Dhandapani, PhD², Philip Leong², and Bryan Tackett, PhD² ¹Weck Laboratories, Inc., 14859 Clark Ave, Hacienda Heights, CA 91745 USA ²Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Overview

EPA method 521 is a popular testing method for nitrosamine compounds by GC-MS/MS. Here, a method that employs activated carbon for the extraction of polar analytes such as N-Nitrosodimethylamine is used. While dealing with polar molecules like nitrosamines, traditional sample preparation using Solid Phase Extraction (SPE) might not be efficient. Since the analyte is polar, a specialized extraction technique is needed. Activated carbon contains porous carbon with defined pore volume and high surface area that is an appropriate choice for polar analytes. In this application note, we have extracted several nitrosamines from water using Strata Activated Carbon cartridges. In addition, multiple batches of Strata Activated Carbon cartridges were evaluated to prove the consistency of the newly developed porous activated carbon material.

Sample Preparation

Condition:	Strata Activated Carbon, 2 g/6 mL cartridge (Part No. <u>CSO-9209</u>) with 2 washes of 3 mL Methylene Chloride, soak for 2 min, 2 washes of 3 mL Methanol, soak 2 min	
Equilibrate:	Cartridges with 3 washes of 3 mL Water	
Load:	500 mL water sample spiked with internal standard onto cartridges	
Dry:	Cartridges for 10 min	
Elute:	With 3 mL Methylene Chloride and repeat elution 3 more times (12 mL total volume)	
Water Removal:	Pass eluent through Sodium Sulfate Giga [™] tubes, 5 g/20 mL (Part No. 8B-5124-LEG), pre-wetted with Methylene Chloride and wash with 3 mL Methylene Chloride	
Evaporate:	Solvent to approximately 0.9 mL in a water bath at 20-25 °C under a gentle stream of Nitrogen	
Adjust:	Final volume to 1 mL with Methylene Chloride	

GC-MS Conditions

Part No.: 7HG-G030-22

Injection Program: Ramp(°C/min) Temp (°C) Time(min)

 40
 0.10

 200
 250
 17.0

Recommended Liner: Zebron PLUS Z-Liner™ (Compatible with Agilent® & Thermo

Scientific® GC instrument)

Liner Part No.: AG2-0A03-05

Carrier Gas: Helium @ 1 mL/min (Constant Flow)

Oven Program: Ramp(°C/min) Temp (°C) Time(min)

- 40 2.5 10 105 0.0 10 200 0.0 50 280 1.5

Detector: MSD

Detector Temperature: 300 °C

Table 1. Mass Spec Transitions

Peak No.	Analyte Name	m/z
1	N-Nitrosodimethylamine-d6 (NDMA-d6) (SS)	81, 59
2	N-Nitrosodimethylamine (NDMA)	75, 56
3	N-Nitrosomethylethylamine (NMEA)	90, 61
4	N-Nitrosodiethylamine-d10 (NDEA-d10) (IS)	114, 81
5	N-Nitrosodiethylamine (NDEA)	85, 75
6	N-Nitrosopropylamine-d14 (NDPA-d14) (IS)	146, 97
7	N-Nitrosodipropylamine (NDPA)	132, 89
8	N-Nitrosopyrrolidine (NPYR)	102, 55
9	N-Nitrosomorpholine (NMOR)	87, 86
10	N-Nitrosopiperidine (NPIP)	116, 69
11	N-Nitrosodibutylamine (NDBA)	159, 103

Figure 1. Strata™ Activated Carbon Extraction Cartridges






Figure 2. Nitrosamine Standards in Water after Strata Activated Carbon Extraction




Peak No.	Analyte	RT (min)
1	NDMA	7.299
2	NMEA	8.967
3	NDEA	10.367
4	NDPA	13.503
5	NPYR	13.551
6	NMOR	13.631
7	NPIP	14.407
8	NDBA	16.623

Figure 3. Nitrosamine Internal Standard and Surrogate Standard Spike in Water after Strata™ Activated Carbon Extraction



Peak No.	Analyte Name	RT (min)
1	NDMA-d6 (SS)	7.236
2	NDEA-d10 (IS)	10.268
3	NDPA-d14 (IS)	13.413

Figure 4. Water Spiked with 2 ppb Nitrosamine Standards, Surrogate Standard, and Internal Standard after Strata **Activated Carbon Extraction**



Peak No.	Analyte	RT (min)
1	NDMA-d6 (SS)	7.236
2	NDMA	7.299
3	NMEA	8.967
4	NDEA-d10 (IS)	10.268
5	NDEA	10.367
6	NDPA-d14 (IS)	13.413
7	NDPA	13.503
8	NPYR	13.551
9	NMOR	13.631
10	NPIP	14.407
11	NDBA	16.623

Figure 5. % Recovery of Surrogate Standard After Strata™ Activated Carbon Extraction




Figure 6. % Recovery of Nitrosamines from a Water Sample After Strata Activated Carbon Extraction





Figure 7. Nitrosamine Recovery Variability After Strata Activated Carbon Extraction



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia t: +61 (0)2-9428-6444 auinfo@phenomenex.com

Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681 info@phenomenex.com

China

t: +86 400-606-8099 cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077 cz-info@phenomenex.com

t: +45 4824 8048 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

Germany t: +49 (0)6021-58830-0 anfrage@phenomenex.com

Hong Kong t: +852 6012 8162 hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400 indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707 indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405 eireinfo@phenomenex.com

Italy t: +39 051 6327511 italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262 ipinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700 nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226 tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700 nlinfo@phenomenex.com

New Zealand t: +64 (0)9-4780951 nzinfo@phenomenex.com

t: +47 810 02 005 nordicinfo@phenomenex.com

t: +48 22 104 21 72 pl-info@phenomenex.com

Portugal

t: +351 221 450 488 ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944 sginfo@phenomenex.com

Slovakia t: +420 272 017 077 sk-info@phenomenex.com

Spain t: +34 91-413-8613 espinfo@phenomenex.com

t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287 thaiinfo@phenomenex.com

United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555 info@phenomenex.com

All other countries/regions Corporate Office USA

t: +1 (310) 212-0555 info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions

Trademarks

Strata, Zebron, 5MS Rus, Z-Liner, Guardian, Giga, and BE-HAPPY are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. Thermo Scientific is a registered of Thermo Fisher Scientific, Inc.

Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Agilent Technologies, Inc. or Thermo Fisher Scientific, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures. © 2022 Phenomenex, Inc. All rights reserved.

