




**KINETEX**<sup>®</sup>  
Core-Shell Technology

# NEW Kinetex **F5**

**HPLC/UHPLC Core-Shell Columns**

Whoa! I can  
even separate  
structural  
isomers!



- Reduce Method Development Time by Days
- Greater Reproducibility than other **PFPs**
- 5 Glorious Interaction Mechanisms
- 5 Valuable LC Separation Modes

F5

TM


**phenomenex**<sup>®</sup>  
...breaking with tradition<sup>SM</sup>



# How I Work

With the astonishing combination of core-shell performance and 5 interaction mechanisms, Kinetex® F5 columns will effortlessly drive your orthogonal HPLC/UHPLC development!

- A Hydrophobic**  
Carbon skeleton of linker and ring encourage neutral/hydrophobic retention
- B Aromatic**  
In non-acetonitrile mobile phases,  $\pi$ - $\pi$  electrons of the carbon ring interact with analyte  $\pi$ - $\pi$  electrons and result in positive retention increase



- C Electrostatic**  
High electronegativity of the fluorine groups create dipole moments, aiding in polar compound retention. Induced dipole moments can also aid neutral compound retention.
- D Steric/Planar**  
Shape selectivity allows for isomeric separations that are otherwise impossible on traditional alkyl phases.
- E Hydrogen Bonding**  
Extremely effective retention mechanism caused as polar functional groups of analyte interact with the electron greedy fluorine.

© 2018 Phenomenex, Inc. All rights reserved.

# Why I'm Better

While older pentafluorophenyl phases (PFP, PFPP, F5, etc.) are based on existing bonding techniques and technologies that promote irreproducibility, the Kinetex F5 was meticulously designed by Phenomenex R&D and its customers, to provide consistently accurate and high performance results.


Learn how the new Kinetex F5 will get you the results you deserve time and time again!



- Core-Shell Advantage ..... **p. 4**
- Method Development Versatility ..... **p. 5**
- Unmatched Reproducibility ..... **pp. 6-7**
- 100% Aqueous Stability ..... **p. 8**
- Isomeric Separations ..... **p. 9**
- Vitamin D3 Epimers ..... **p. 10**
- Multi Component Screening ..... **p. 11**
- Methanol vs. Acetonitrile ..... **p. 12**
- Trace Impurity Detection ..... **p. 13**
- UHPLC Consistency  
and Performance ..... **pp. 14-15**
- Kinetex Column Selection ..... **pp. 16-17**
- Ordering Information ..... **pp. 18-21**



If you are not completely satisfied with Kinetex core-shell columns, send in your comparative data to a similar product with the Kinetex column within 45 days for a FULL REFUND.




# Why Wait?


## You Can Get Improvements Today!

The band broadening (wide peaks) and lengthy retention times of traditional fully porous products can be limiting your results. Now you can use the incredible efficiency levels of Kinetex® Core-Shell Technology to achieve shorter run times, higher levels of sensitivity, and overall better HPLC/UHPLC results.

### Core-Shell vs. Fully Porous Efficiency Levels (plates/m)




### Core-Shell Performance Gains



# Method Development

## Versatility

Combine core-shell performance, multiple retention mechanisms and the Kinetex F5 column's ability to be run in a variety of separation modes and you now have an impeccable method development tool at your disposal.






# Dependability




## Batch-to-Batch, Column-to-Column

Conventional fully porous and core-shell PFP/F5 columns fail to reach the level of repeatability that you deserve. Inconsistencies in their base silica have led to data inaccuracies that waste your time and money. Kinetex® F5 columns were specifically designed to avoid these past problems and provide you with the highest degree of reproducibility.

### Kinetex 2.6µm F5



### Supelco® Ascentis® Express 2.7µm F5



**Conditions for all columns:**

**Column:** Kinetex 2.6µm F5  
 Discovery 3µm HS F5  
 Ascentis Express 2.7µm F5  
 Hypersil GOLD 3µm PFP  
 Accucore 2.6µm PFP

**Dimensions:** 50 x 4.6 mm

**Mobile Phase:** A: 0.1 % Formic acid in Water  
 B: 0.1 % Formic acid in Acetonitrile

**Gradient:** 5-95 % B over 5 minutes

**Flow Rate:** 1.85 mL/min

**Temperature:** Ambient


**Detection:** UV @ 254 nm

- Sample:**
1. Uracil
  2. Pindolol
  3. Chlorpheniramine
  4. Nortriptyline
  5. 3-Methyl-4-Nitrobenzoic acid
  6. 5-Methyl Salicyl Aldehyde
  7. Hexaphenone


What happened to peaks 3 and 4?

Accucore is a trademark and Hypersil GOLD is registered trademark of Thermo Hypersil-Keystone. Ascentis, Discovery and Supelco are registered trademarks of Sigma-Aldrich Co. LLC. Phenomenex is not affiliated with any of the above companies. Comparative separations may not be representative of all applications.

### Supelco® Discovery® 3µm HS F5




### Thermo Hypersil GOLD® 3µm PFP



**Adsorption of Peaks 3 and 4!**

**Retention Time Shifting**

### Thermo Accucore™ 2.6µm PFP



**Conditions for all columns:**

**Column:** Kinetex 2.6µm F5  
Discovery 3µm HS F5  
Ascentis Express 2.7µm F5  
Hypersil GOLD 3µm PFP  
Accucore 2.6µm PFP

**Dimensions:** 50 x 4.6 mm

**Mobile Phase:** A: 0.1% Formic acid in Water  
B: 0.1% Formic acid in Acetonitrile

**Gradient:** 5-95% B over 5 minutes

**Flow Rate:** 1.85 mL/min


**Temperature:** Ambient

**Detection:** UV @ 254 nm

**Sample:** 1. Uracil  
2. Pindolol  
3. Chlorpheniramine  
4. Nortriptyline  
5. 3-Methyl-4-Nitrobenzoic acid  
6. 5-Methyl Salicyl Aldehyde  
7. Hexaphenone

**Are you okay with peak adsorptions and retention time shifts?**




Accucore is a trademark and Hypersil GOLD is registered trademark of Thermo Hypersil-Keystone. Ascentis, Discovery and Supelco are registered trademarks of Sigma-Aldrich Co. LLC. Phenomenex is not affiliated with any of the above companies. Comparative separations may not be representative of all applications.



# Extra Polar Retention and 100% Aqueous Stability



Unlike traditional alkyl stationary phases, the polar functional group of the Kinetex® F5 makes it stable in 100% aqueous mobile phase conditions. This can be especially beneficial for methods that require the retention and resolution of problematic polar compounds that typically have poor retention under reversed phase conditions.

## Organic Acids



Each successive day, column was equilibrated, run and stored in the 100% aqueous conditions mentioned above.

**Column:** Kinetex 2.6µm F5  
**Dimensions:** 100 x 4.6mm  
**Part No.:** 00D-4723-E0  
**Mobile Phase:** 20 mM Sodium phosphate pH 2.5  
**Flow Rate:** 1.5 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 210 nm  
**Sample:** 1. Tartaric Acid  
2. Malic Acid  
3. Maleic Acid  
4. Fumaric Acid  
5. Succinic Acid



**INCREDIBLE!**  
Even the isomers  
Maleic acid and  
Fumaric acid  
are separated!






# Novel Selectivity For Isomeric Separations

While a C18 can differentiate between the small addition of a single methyl group between two similar compounds, it cannot separate compounds with only structural differences, like positional isomers. This is where the electrostatic and planar interactions of the Kinetex F5 give focused attention to the resolution you need.

## Methoxybenzene Isomers



Conditions for all columns:

**Column:** Kinetex 2.6µm F5  
 Kinetex 2.6µm C18  
 Kinetex 2.6µm Biphenyl

**Dimensions:** 150 x 4.6 mm

**Mobile Phase:** A: 0.1% TFA in Water  
 B: Acetonitrile



**Isocratic:** A/B (65:35)

**Flow Rate:** 1 mL/min

**Temperature:** Ambient

**Detection:** UV @ 254 nm

**Sample:** 1. 1,2,3-Trimethoxybenzene  
 2. 1,2-Dimethoxybenzene  
 3. 1,2,4-Trimethoxybenzene  
 4. 1,4-Dimethoxybenzene  
 5. Methoxybenzene  
 6. 1,3-Dimethoxybenzene




# Clinical Research

## Vitamin D3 Epimers

Even tandem mass spectrometry (LC/MS/MS) analysis has need to utilize the wonderful cross functionality of Kinetex® F5. With the same fragment ions coming from the Vitamin D3 epimers, reproducible chromatographic separation is a must. The unique combination of polar/non-polar selectivity of Kinetex F5 easily performs this necessary separation in a highly sensitive and short analysis window.


### Kinetex 2.6 µm F5




#### Conditions for all columns:

- Column:** Kinetex 2.6 µm F5  
HALO 2.7 µm PFP  
XSelect HSS 2.5 µm PFP
- Dimensions:** 100 x 4.6 mm
- Mobile Phase:** A: 0.1% Formic acid in Water  
B: 0.1% Formic acid in Methanol
- Isocratic:** A/B (15:85)
- Flow Rate:** 0.75 mL/min
- Temperature:** Ambient
- Detection:** MS/MS (SCIEX API 4000™)
- Sample:** 1. 25-OH Vitamin D3  
2. 25-OH Vitamin D2  
3. 3-epi-25-OH Vitamin D3

### Advanced Materials Technology HALO® 2.7 µm PFP



### Waters® XSelect® HSS 2.5 µm PFP



Vitamin D3 epimers are not fully separated


HALO is a registered trademark of Advanced Materials Technology, Inc. Waters and XSelect are registered trademarks of Waters Corporation. Phenomenex is not affiliated with any of the above companies. Comparative separations may not be representative of all applications.

# Food Testing

# Multi Component Analysis


The versatility of Kinetex F5 core-shell columns matches marvelously with ingredient quantification and contaminant identification. Acidic food additives may tail and lack retention on a C18 column, but the Kinetex F5 offers superb polar retention and performance. Meanwhile, large multi component contaminant screens and their mixtures of acids, bases, neutrals and structurally similar compounds, can easily utilize the combination of polar, non-polar and geometric interaction mechanisms of the Kinetex F5 to get excellent separation and sensitivity.

## Food Additives



**Column:** Kinetex 2.6µm F5  
**Dimensions:** 150 x 2.1 mm  
**Part No.:** 00F-4723-AN  
**Mobile Phase:** A: 0.1% Phosphoric acid in Water  
 B: Acetonitrile  
**Gradient:** 5-35% A in 6 min. Hold for 4 min.  
**Flow Rate:** 0.6 mL/min  
**Temperature:** 30 °C  
**Detection:** UV @ 240 nm  
**Sample:** 1. Saccharin  
 2. p-Hydroxybenzoic acid  
 3. Sorbic acid  
 4. Dehydroacetic acid  
 5. Methyl paraben

## Antibiotics Screen



**Column:** Kinetex 2.6µm F5  
**Dimensions:** 50 x 2.1 mm  
**Part No.:** 00B-4723-AN  
**Mobile Phase:** A: 0.1% Formic acid in Water  
 B: 0.1% Formic acid in Methanol  

Gradient:	Time (min)	% B
	0	0
	0.08	0
	1.08	100
	8.08	100
	11.08	0
	11.58	0

**Flow Rate:** 0.5 mL/min  
**Temperature:** 30 °C  
**Detection:** MS/MS (SCIEX API 4000™)  
**Sample:** 1. Amoxicillin  
 2. Cefalexin  
 3. Cefazolin  
 4. Cefoperazone  
 5. Cefapirin  
 6. Cloxacillin  
 7. Dicloxacillin  
 8. Ciprofloxacin  
 9. Difloxacin  
 10. Marbofloxacin  
 11. Sulfadiazine  
 12. Sulfamerazine  
 13. Sulfamethazine  
 14. Sulfamethoxazole  
 15. Sulfapyridine  
 16. Sulfaquinoxaline  
 17. Sulfathiazole  
 18. Neospiramycin  
 19. Doxycycline  
 20. Tiamulin  
 21. Valnemulin  
 22. Rifaximin  
 23. Lincomycin  
 24. Nafcillin

Find more Kinetex F5 applications at [www.phenomenex.com/KinetexF5Apps](http://www.phenomenex.com/KinetexF5Apps)





# Methanol vs. Acetonitrile

## A Phenyl Story

While mobile phase modifiers can help adjust retention, with the use of the Kinetex® F5, the major mobile phase constituent can also be used to manipulate elution order and retention properties. Acetonitrile can be used to disrupt pi-pi interactions between compounds and phenyl phases, while switching to the weaker solvent Methanol will encourage aromatic interactions.

### Methoxybenzene Isomers



#### Conditions for all columns:

**Column:** Kinetex 2.6  $\mu$ m F5

**Dimensions:** 150 x 4.6 mm

**Part No.:** 00F-4723-E0

**Mobile Phase:** as noted



**Isocratic:** A/B (65:35)

**Flow Rate:** 1 mL/min

**Temperature:** Ambient

**Detection:** UV @ 254 nm

**Sample:** 1. 1,2,3-Trimethoxybenzene  
2. 1,2-Dimethoxybenzene  
3. 1,2,4-Trimethoxybenzene




**Impact selectivity and retention with solvent adjustment!**


# Trace Pharmaceutical Impurity Detection

Trace impurities of active pharmaceutical ingredients are incredibly important to identify and quantify. With the rapid performance value of core-shell technology combined with the versatility of a pentafluorophenyl, the Kinetex F5 is the precise alternative to other reversed phase columns that you need. Easily utilize the Kinetex F5 to get greater sensitivity, better resolution and all in shorter analysis times.

**Kinetex 1.7  $\mu$ m F5**



**Waters® ACQUITY® CSH™ 1.7  $\mu$ m Fluoro-Phenyl**



**Conditions for all columns:**

**Column:** Kinetex 1.7  $\mu$ m F5  
 ACQUITY CSH 1.7  $\mu$ m Fluoro-Phenyl

**Dimensions:** 50 x 2.1 mm

**Mobile Phase:** A: 20mM Potassium phosphate pH 2.3  
 B: Methanol

Gradient:	Time (min)	% B
	0	5
	10	95
	10.01	5

**Flow Rate:** 0.3 mL/min

**Temperature:** Ambient


**Detection:** UV @ 254 nm

**Sample:**  
 1. Impurity 1  
 2. Impurity 2  
 3. Proprietary Active Pharmaceutical Ingredient



# Incredible UHPLC Reproducibility

Highly sensitive UHPLC separations are dependent upon consistent quantitation and consistent results. With the reliability of the Kinetex® F5, you no longer have to settle for the inconsistent results of PFP products that currently exist on the market. Our highest standards of quality will ensure that you are fully satisfied with each and every Kinetex F5 column!

## Trusted Consistency



**Column:** Kinetex 1.7  $\mu$ m F5  
**Dimensions:** 50 x 4.6 mm  
**Mobile Phase:** A: 0.1 % Formic in Water  
B: 0.1 % Formic in Acetonitrile  
**Gradient:** 5-95 % B over 5 min.  
**Flow Rate:** 1.85 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 254 nm  
**Sample:** 1. Uracil  
2. Pindolol  
3. Chlorpheniramine  
4. Nortriptyline  
5. 3-Methyl-4-Nitrobenzoic acid  
6. 5-Methyl Salicyl Aldehyde  
7. Hexaphenone




# A Better PFP

## Selectivity for UHPLC

Kinetex 1.7  $\mu\text{m}$  core-shell technology produces increased efficiencies over traditional sub-2  $\mu\text{m}$  columns on the market, yielding remarkable chromatographic resolution, higher peak capacities, and greater sensitivity. Add these performance gains alongside the novel selectivity and excellent reproducibility of the Kinetex F5 and you now have an incredible UHPLC solution at your fingertips.

### Pharmaceutical Mix



**Conditions for all columns:**

**Column:** Kinetex 1.7  $\mu\text{m}$  F5  
ACQUITY CSH 1.7  $\mu\text{m}$  Fluoro-Phenyl  
ACQUITY HSS 1.8  $\mu\text{m}$  PFP

**Dimensions:** 50 x 2.1 mm

**Mobile Phase:** A: 10 mM Ammonium acetate pH 3.2  
B: Acetonitrile



Gradient:	Time (min)	% B
	90	10
	90	10
	5	95
	5	95
	90	10
	90	10

**Flow Rate:** 0.4 mL/min

**Temperature:** Ambient

**Detection:** UV @ 220 nm

**Sample:**  
1. Pindolol  
2. Propranolol  
3. Indoprofen  
4. Naproxen  
5. Warfarin  
6. Terfenidine



Terrellasaurus – first discovered in Torrance, CA.

Waters and ACQUITY are registered trademarks, and CSH is a trademark of Waters Corporation. Phenomenex is not affiliated with Waters Corporation. Comparative separations may not be representative of all applications.

# Get Method Development Flexibility For Your Small Molecule Analysis

## Recommended Selectivities If You're Working With:

### Acids

- Kinetex® C18
- Kinetex F5
- Kinetex Phenyl-Hexyl

### Bases

- Kinetex EVO C18
- Kinetex XB-C18
- Kinetex Biphenyl
- Kinetex Polar C18

### Neutrals

- Kinetex C18
- Kinetex C8
- Kinetex Biphenyl

### Aromatics

- Kinetex Biphenyl
- Kinetex Phenyl-Hexyl
- Kinetex F5

### Acids, Bases, and Neutrals

- Kinetex Polar C18
- Kinetex Biphenyl
- Kinetex EVO C18
- Kinetex F5

### Highly Polar Compounds

- Kinetex Polar C18
- Kinetex F5
- Kinetex Biphenyl
- Kinetex HILIC

### High pH

- Kinetex EVO C18

### Isomers

- Kinetex F5

## Upgrading Your Fully Porous Methods:

### Fully Porous 3 $\mu\text{m}$ or 5 $\mu\text{m}$

- **Kinetex 5  $\mu\text{m}$**  – Drop-in for easy performance improvements with no backpressure increase
- **Kinetex 3.5  $\mu\text{m}$**  – Drop-in for easy performance improvements of pharmacopoeia methods
- **Kinetex 2.6  $\mu\text{m}$**  – Dramatically improve results with efficiency/peak capacity gains

### Fully Porous sub-2 $\mu\text{m}$

- **Kinetex 2.6  $\mu\text{m}$**  – Get similar efficiencies at lower backpressure allowing for greater productivity gains
- **Kinetex 1.7  $\mu\text{m}$**  – Up to 20% greater efficiencies resulting in drop-in improvements
- **Kinetex 1.3  $\mu\text{m}$**  – Incredible efficiency gains on high end UHPLC systems

### Fully Porous Preparative LC

- **Kinetex 5  $\mu\text{m}$**  – Drop-in for easy performance improvement with no backpressure increase



# Choosing The Best Core-Shell Platform For You is Easy!

## For Small Molecules

	5 $\mu$ m	3.5 $\mu$ m	2.6 $\mu$ m	1.7 $\mu$ m	1.3 $\mu$ m
<b>UHPLC</b>					
<b>HPLC</b>					
<b>PREP LC</b>					

Phase	Best Use	pH Stability	Available Particle Size(s)				
<b>F5</b>	Highly reproducible pentafluorophenyl propyl phase that offers a unique combination of polar, hydrophobic, aromatic, and shape selectivity	1.5 - 8.5*	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>Polar C18</b>	C18 provides all purpose non-polar interactions, while novel polar modified surface increases polar compound retention and provides 100% aqueous stability	1.5 - 8.5*			2.6 $\mu$ m		
<b>EVO C18</b>	Robust reversed phase methods even in alkaline conditions with improved peak shape for polar basic compounds	1 - 12	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>C18</b>	All purpose phase that offers the hydrophobic retention and methylene selectivity chromatographers expect from a C18 column	1.5 - 8.5*	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	1.3 $\mu$ m
<b>XB-C18</b>	C18 phase with protective butyl side chains for improved peak shape for basic compounds under neutral and acidic conditions	1.5 - 8.5*	5 $\mu$ m	3.5 $\mu$ m	2.6 $\mu$ m	1.7 $\mu$ m	
<b>C8</b>	USP L7 phase that provides less hydrophobic and methylene selectivity than a C18	1.5 - 8.5*	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>Biphenyl</b>	100% aqueous stable and allows for excellent reversed phase retention and enhanced polar and aromatic selectivity	1.5 - 8.5*	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>Phenyl-Hexyl</b>	Reversed phase chemistry that allows for greater retention and separation of aromatic hydrocarbons	1.5 - 8.5*	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>HILIC</b>	Unbonded silica phase for HILIC conditions to provide selectivity for polar compounds	2.0 - 7.5	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	

\*pH stability under gradient conditions. pH stability is 1.5-10 under isocratic conditions.


## Phenomenex Application Specific Core-Shell Products

Material	Phase	Best Use	pH Stability	Available Particle Size(s)				
<b>For Peptides (<math>\leq</math> 10,000 Da)</b>								
<b>Aeris™ PEPTIDE</b>	XB-C18	Excellent hydrophobicity and methylene selectivity for peptide and peptide mapping separations	1.5 - 9.0	5 $\mu$ m	3.6 $\mu$ m	2.6 $\mu$ m	1.7 $\mu$ m	
<b>For Proteins (&gt; 10,000 Da)</b>								
<b>Aeris WIDEPORÉ</b>	XB-C18	Maximum hydrophobicity and high temp stability for hydrophilic and PEGylated proteins	1.5 - 9.0		3.6 $\mu$ m			
	XB-C8	Medium hydrophobicity and high temp stability for moderately hydrophobic proteins and glycosylated proteins	1.5 - 9.0		3.6 $\mu$ m			
	C4	Lowest hydrophobicity for very large or very hydrophobic proteins	1.5 - 9.0		3.6 $\mu$ m			
<b>For Synthetic Oligonucleotides (DNA/RNA)</b>								
<b>Clarity® Oligo-XT</b>	C18	Rapid, high efficiency reversed phase LC/MS analysis for QC and characterization	1 - 12	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	

# Protect Any UHPLC Column

Protect your UHPLC column, including Kinetex core-shell columns, from damaging contaminants and microparticulates with the SecurityGuard ULTRA guard cartridge system!

- Simple to use
- Extend column lifetime
- Pressure rated to 20,000 psi (1,378 bar)
- Fits virtually all manufacturers' columns 2.1 to 4.6 mm ID



See it in action:

[www.phenomenex.com/SecurityGuardULTRA](http://www.phenomenex.com/SecurityGuardULTRA)

# Ordering Information

5 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges <sup>‡</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>EVO C18</b>	00A-4633-AN	00B-4633-AN	00D-4633-AN	00F-4633-AN	AJ0-9298
<b>F5</b>	00A-4724-AN	00B-4724-AN	00D-4724-AN	00F-4724-AN	AJ0-9322
<b>Biphenyl</b>	00A-4627-AN	00B-4627-AN	00D-4627-AN	—	AJ0-9209
<b>XB-C18</b>	00A-4605-AN	00B-4605-AN	00D-4605-AN	—	AJ0-8782
<b>C18</b>	00A-4601-AN	00B-4601-AN	00D-4601-AN	00F-4601-AN	AJ0-8782
<b>C8</b>	—	00B-4608-AN	00D-4608-AN	—	AJ0-8784
<b>Phenyl-Hexyl</b>	—	00B-4603-AN	00D-4603-AN	—	AJ0-8788

for 2.1 mm ID

5 µm MidBore™ Columns (mm)				SecurityGuard ULTRA Cartridges <sup>‡</sup>
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
<b>EVO C18</b>	00B-4633-Y0	00D-4633-Y0	00F-4633-Y0	AJ0-9297
<b>F5</b>	00B-4724-Y0	00D-4724-Y0	00F-4724-Y0	AJ0-9321
<b>Biphenyl</b>	00B-4627-Y0	00D-4627-Y0	00F-4627-Y0	AJ0-9208
<b>XB-C18</b>	00B-4605-Y0	00D-4605-Y0	00F-4605-Y0	AJ0-8775
<b>C18</b>	00B-4601-Y0	00D-4601-Y0	00F-4601-Y0	AJ0-8775
<b>C8</b>	00B-4608-Y0	00D-4608-Y0	—	AJ0-8777
<b>Phenyl-Hexyl</b>	00B-4603-Y0	00D-4603-Y0	—	AJ0-8781

for 3.0 mm ID

5 µm Analytical Columns (mm)					SecurityGuard ULTRA Cartridges <sup>‡</sup>
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
<b>EVO C18</b>	00B-4633-E0	00D-4633-E0	00F-4633-E0	00G-4633-E0	AJ0-9296
<b>F5</b>	00B-4724-E0	00D-4724-E0	00F-4724-E0	00G-4724-E0	AJ0-9320
<b>Biphenyl</b>	00B-4627-E0	00D-4627-E0	00F-4627-E0	00G-4627-E0	AJ0-9207
<b>XB-C18</b>	00B-4605-E0	00D-4605-E0	00F-4605-E0	00G-4605-E0	AJ0-8768
<b>C18</b>	00B-4601-E0	00D-4601-E0	00F-4601-E0	00G-4601-E0	AJ0-8768
<b>C8</b>	00B-4608-E0	00D-4608-E0	00F-4608-E0	00G-4608-E0	AJ0-8770
<b>Phenyl-Hexyl</b>	00B-4603-E0	00D-4603-E0	00F-4603-E0	00G-4603-E0	AJ0-8774

for 4.6 mm ID

5 µm Semi-Preparative Columns (mm)			SecurityGuard SemiPrep Cartridges <sup>***</sup>
Phases	150 x 10	250 x 10	10 x 10
<b>EVO C18</b>	00F-4633-N0	00G-4633-N0	AJ0-9306 /3pk
<b>F5</b>	—	00G-4724-N0	AJ0-9323 /3pk
<b>C18</b>	00F-4601-N0	00G-4601-N0	AJ0-9278
<b>Biphenyl</b>	00F-4627-N0	00G-4627-N0	AJ0-9280

for 10 mm ID

5 µm Axia™ Packed Preparative Columns (mm)					SecurityGuard PREP Cartridges <sup>*</sup>
Phases	50 x 21.2	100 x 21.2	150 x 21.2	250 x 21.2	15 x 21.2
<b>EVO C18</b>	00B-4633-P0-AX	00D-4633-P0-AX	00F-4633-P0-AX	00G-4633-P0-AX	AJ0-9304 /ea
<b>F5</b>	—	—	00F-4724-P0-AX	00G-4724-P0-AX	AJ0-9324
<b>Biphenyl</b>	00B-4627-P0-AX	00D-4627-P0-AX	00F-4627-P0-AX	00G-4627-P0-AX	AJ0-9272
<b>XB-C18</b>	00B-4605-P0-AX	00D-4605-P0-AX	00F-4605-P0-AX	00G-4605-P0-AX	AJ0-9145
<b>C18</b>	00B-4601-P0-AX	00D-4601-P0-AX	00F-4601-P0-AX	00G-4601-P0-AX	AJ0-9145
<b>C8</b>	00B-4608-P0-AX	00D-4608-P0-AX	00F-4608-P0-AX	00G-4608-P0-AX	AJ0-9205
<b>Phenyl-Hexyl</b>	00B-4603-P0-AX	00D-4603-P0-AX	00F-4603-P0-AX	00G-4603-P0-AX	AJ0-9147
<b>HILIC</b>	—	00D-4606-P0-AX	00F-4606-P0-AX	00G-4606-P0-AX	AJ0-9277

for 21.2 mm ID

<sup>‡</sup> SecurityGuard ULTRA Cartridges require holder, Part No.: AJ0-9000

<sup>\*</sup> PREP SecurityGuard Cartridges require holder, Part No.: AJ0-8223

<sup>\*\*</sup> PREP SecurityGuard Cartridges require holder, Part No.: AJ0-8277

<sup>\*\*\*</sup> SemiPrep SecurityGuard Cartridges require holder, Part No.: AJ0-9281

# Ordering Information

5 µm Axia Packed Preparative Columns (mm)					SecurityGuard PREP Cartridges**
Phases	50 x 30	100 x 30	150 x 30	250 x 30	15 x 30
					/ea
<b>EVO C18</b>	00B-4633-U0-AX	00D-4633-U0-AX	00F-4633-U0-AX	00G-4633-U0-AX	AJO-9305
					/ea
<b>F5</b>	00B-4724-U0-AX	00D-4724-U0-AX	00F-4724-U0-AX	00G-4724-U0-AX	AJO-9325
<b>Biphenyl</b>	—	—	00F-4627-U0-AX	—	AJO-9273
<b>XB-C18</b>	00B-4605-U0-AX	00D-4605-U0-AX	00F-4605-U0-AX	00G-4605-U0-AX	AJO-9204
<b>C18</b>	00B-4601-U0-AX	00D-4601-U0-AX	00F-4601-U0-AX	00G-4601-U0-AX	AJO-9204
<b>C8</b>	00B-4608-U0-AX	00D-4608-U0-AX	00F-4608-U0-AX	00G-4608-U0-AX	AJO-9217
<b>Phenyl-Hexyl</b>	00B-4603-U0-AX	00D-4603-U0-AX	00F-4603-U0-AX	00G-4603-U0-AX	AJO-9216

for 30 mm ID

3.5 µm Analytical Columns (mm)			SecurityGuard ULTRA Cartridges†
Phases	100 x 4.6	150 x 4.6	3/pk
<b>XB-C18</b>	00D-4744-E0	00F-4744-E0	AJO-8768

for 4.6 mm ID

2.6 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
<b>XB-C18</b>	00B-4496-A0	00D-4496-A0	00F-4496-A0

2.6 µm Minibore Columns (mm)						SecurityGuard™ ULTRA Cartridges†
Phases	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>EVO C18</b>	00A-4725-AN	00B-4725-AN	—	00D-4725-AN	00F-4725-AN	AJO-9298
<b>Polar C18</b>	00A-4759-AN	00B-4759-AN	—	00D-4759-AN	00F-4759-AN	AJO-9532
<b>F5</b>	00A-4723-AN	00B-4723-AN	—	00D-4723-AN	00F-4723-AN	AJO-9322
<b>Biphenyl</b>	00A-4622-AN	00B-4622-AN	—	00D-4622-AN	00F-4622-AN	AJO-9209
<b>XB-C18</b>	00A-4496-AN	00B-4496-AN	00C-4496-AN	00D-4496-AN	00F-4496-AN	AJO-8782
<b>C18</b>	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJO-8782
<b>C8</b>	00A-4497-AN	00B-4497-AN	00C-4497-AN	00D-4497-AN	00F-4497-AN	AJO-8784
<b>HILIC</b>	00A-4461-AN	00B-4461-AN	00C-4461-AN	00D-4461-AN	00F-4461-AN	AJO-8786
<b>Phenyl-Hexyl</b>	00A-4495-AN	00B-4495-AN	00C-4495-AN	00D-4495-AN	00F-4495-AN	AJO-8788

for 2.1 mm ID

2.6 µm MidBore™ Columns (mm)						SecurityGuard ULTRA Cartridges†
Phases	30 x 3.0	50 x 3.0	75 x 3.0	100 x 3.0	150 x 3.0	3/pk
<b>EVO C18</b>	—	00B-4725-Y0	—	00D-4725-Y0	00F-4725-Y0	AJO-9297
<b>Polar C18</b>	—	00B-4759-Y0	—	00D-4759-Y0	00F-4759-Y0	AJO-9531
<b>F5</b>	—	00B-4723-Y0	—	00D-4723-Y0	00F-4723-Y0	AJO-9321
<b>Biphenyl</b>	—	00B-4622-Y0	—	00D-4622-Y0	00F-4622-Y0	AJO-9208
<b>XB-C18</b>	00A-4496-Y0	00B-4496-Y0	00C-4496-Y0	00D-4496-Y0	00F-4496-Y0	AJO-8775
<b>C18</b>	00A-4462-Y0	00B-4462-Y0	00C-4462-Y0	00D-4462-Y0	00F-4462-Y0	AJO-8775
<b>C8</b>	00A-4497-Y0	00B-4497-Y0	00C-4497-Y0	00D-4497-Y0	00F-4497-Y0	AJO-8777
<b>HILIC</b>	00A-4461-Y0	—	—	—	00F-4461-Y0	AJO-8779
<b>Phenyl-Hexyl</b>	—	00B-4495-Y0	—	00D-4495-Y0	00F-4495-Y0	AJO-8781

for 3.0 mm ID

† SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

\* PREP SecurityGuard Cartridges require holder, Part No.: AJO-8223

\*\* PREP SecurityGuard Cartridges require holder, Part No.: AJO-8277

\*\*\* SemiPrep SecurityGuard Cartridges require holder, Part No.: AJO-9281

# Ordering Information

2.6 µm Analytical Columns (mm)						SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	3/pk
<b>EVO C18</b>	—	00B-4725-E0	—	00D-4725-E0	00F-4725-E0	AJO-9296
<b>Polar C18</b>	—	00B-4759-E0	—	00D-4759-E0	00F-4759-E0	AJO-9530
<b>F5</b>	—	00B-4723-E0	—	00D-4723-E0	00F-4723-E0	AJO-9320
<b>Biphenyl</b>	—	00B-4622-E0	—	00D-4622-E0	00F-4622-E0	AJO-9207
<b>XB-C18</b>	—	00B-4496-E0	00C-4496-E0	00D-4496-E0	00F-4496-E0	AJO-8768
<b>C18</b>	00A-4462-E0	00B-4462-E0	00C-4462-E0	00D-4462-E0	00F-4462-E0	AJO-8768
<b>C8</b>	—	00B-4497-E0	00C-4497-E0	00D-4497-E0	00F-4497-E0	AJO-8770
<b>HILIC</b>	—	00B-4461-E0	00C-4461-E0	00D-4461-E0	00F-4461-E0	AJO-8772
<b>Phenyl-Hexyl</b>	—	00B-4495-E0	00C-4495-E0	00D-4495-E0	00F-4495-E0	AJO-8774

for 4.6 mm ID

1.7 µm Minibore Columns (mm)					SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>EVO C18</b>	—	00B-4726-AN	00D-4726-AN	00F-4726-AN	AJO-9298
<b>F5</b>	—	00B-4722-AN	00D-4722-AN	00F-4722-AN	AJO-9322
<b>Biphenyl</b>	—	00B-4628-AN	00D-4628-AN	00F-4628-AN	AJO-9209
<b>XB-C18</b>	00A-4498-AN	00B-4498-AN	00D-4498-AN	00F-4498-AN	AJO-8782
<b>C18</b>	00A-4475-AN	00B-4475-AN	00D-4475-AN	00F-4475-AN	AJO-8782
<b>C8</b>	00A-4499-AN	00B-4499-AN	00D-4499-AN	00F-4499-AN	AJO-8784
<b>HILIC</b>	00A-4474-AN	00B-4474-AN	00D-4474-AN	—	AJO-8786
<b>Phenyl-Hexyl</b>	—	00B-4500-AN	00D-4500-AN	00F-4500-AN	AJO-8788

for 2.1 mm ID

1.7 µm MidBore™ Columns (mm)				SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 3.0	50 x 3.0	100 x 3.0	3/pk
<b>XB-C18</b>	00A-4498-Y0	00B-4498-Y0	00D-4498-Y0	AJO-8775
<b>C18</b>	—	00B-4475-Y0	00D-4475-Y0	AJO-8775
<b>C8</b>	00A-4499-Y0	00B-4499-Y0	00D-4499-Y0	AJO-8777
<b>HILIC</b>	—	00B-4474-Y0	—	AJO-8779

for 3.0 mm ID

1.3 µm Minibore Columns (mm)		
Phases	30 x 2.1	50 x 2.1
<b>C18</b>	00A-4515-AN	00B-4515-AN

<sup>†</sup> SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000.

**Protect Your  
LC Investment!**

**novum™**  
simplified liquid extraction



Remove interferences  
such as proteins,  
phospholipids, salts,  
and more in a quick and  
simple procedure.

[www.phenomenex.com/  
Novum](http://www.phenomenex.com/Novum)

**PHENEX™**  
SYRINGE FILTERS



Convenient and  
inexpensive way to  
remove particulates from  
HPLC samples.

[www.phenomenex.com/  
Phenex](http://www.phenomenex.com/Phenex)

**guarantee**

If you are not completely satisfied with Kinetex core-shell columns, send in your comparative data to a similar product with the Kinetex column within 45 days for a **FULL REFUND**.

#### Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

#### Trademarks

Phenomenex, Clarity, and Kinetex are registered trademarks of Phenomenex. Aeris, Axia, MidBore, Novum, Phenex, and SecurityGuard are trademarks of Phenomenex. Waters, ACQUITY, and XSelect are registered trademarks, and CSH is a trademark of Waters Corporation. Agilent and Poroshell are registered trademarks of Agilent Technologies, Inc. Ascentis, Discovery, and Supelco are registered trademarks of Sigma-Aldrich Co. LLC. Accucore is a trademark and Hypersil GOLD is registered trademark of Thermo Hypersil-Keystone. API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license. HALO is a registered trademark of Advanced Materials Technology, Inc.

#### Disclaimer

Phenomenex is not affiliated with Advanced Materials Technology, Inc., Agilent Technologies, Sigma-Aldrich Co., Thermo Hypersil-Keystone or Waters Corporation.

Comparative separations may not be representative of all applications.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362.

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

Novum is patent pending.


FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2018 Phenomenex, Inc. All rights reserved.

# NEW Kinetex F5

## HPLC/UHPLC Core-Shell Columns

- Reduce Method Development Time by Days
- Greater Reproducibility than other **PFPs**
- 5 Glorious Interaction Mechanisms
- 5 Valuable LC Separation Modes



### Australia

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
info@phenomenex.com

### China

t: +86 400-606-8099  
cninfo@phenomenex.com

### Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

### India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

### Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

### Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

### Singapore

t: +65 800-852-3944  
sginfo@phenomenex.com

### Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

### Switzerland

t: +41 61 692 20 20  
swissinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
info@phenomenex.com

### All other countries Corporate Office USA

t: +1 (310) 212-0555  
info@phenomenex.com